

SENSO Clean

Genuine innovation in the field of infection prevention: the digital scrub-up trough with sensor control



ALL BENEFITS AT A GLANCE

SPOTLIGHTS

The **SENSO** *Clean* is equipped with warm white LEDs. The light comes on automatically with the water. The beam of light points towards the stream of water to ensure an optimum illumination of the hands.

TOUCHLESS TEMPERATURE CONTROL

The user can regulate the temperature continuously via two sensor fields. A preset mechanical 3-way valve with motor drive (thermostat-controlled) ensures smooth water temperature adjustement. The integrated anti-scaling technology does not allow the water temperature to exceed 42°C.

SOAP AND DISINFECTANT DISPENSER

The dispensers are controlled by an electronically regulated pump. Disinfectants and liquid soap can be dosed without any contact. This reduces the risk of an accidental transfer of germs through contact to the hands that are to be cleaned. Thanks to integrated sensors, the levels are displayed prominently on the screens.

SOFTWARE-GUIDED HAND-WASHING

The SENSO Clean has a 10.1" monitor at each washbasin. When the program is activated, the integrated software begins the display of the 6 steps for an optimum hand-washing result incl. the specified time according to DIN standards. This process is extended by 2 preceding steps for surgical hand disinfection. Apart from the graphic animation, a feedback system accompanies the washing process for the hands.

UV LEDS FOR QUALITY CONTROL

The UV-A LEDs implemented in the housing allow a quality control after hand-washing. This function is only provided for training purposes.







SENSO Clean 1400

SENSO Clean INTERACTIVE SOFTWARE*

HAND-WASHING AND SURGICAL HAND DISINFECTION

Surgical hand disinfection according to EN 12791



1ST STEP Rub in disinfectant in individ-

ual portions. Begin by disin-

fecting the hands and lower

arms up to the elbows



1ST STEP Rub palm on palm



2ND STEP Then disinfect half the lower arm and hands

4TH STEP Rub outside of fingers on opposite hands with interlocked fingers

FEEDBACK SYSTEM

A series of graphic illustrations shows the steps specified for hand-washing and/or hand disinfection. An accompanying feedback system shows the individual progress and helps motivate the doctors and clinic staff.



TIMER



A count-down begins automatically at the start of the washing process and the disinfection. This makes it easier for the user to observe the recommended washing time / prescribed disinfection.



COLD

Software guided, fully automated hand disinfection

Each washbasin has a 10.1" monitor and the integrated guiding software serves quality compliance and control

Visual feedback for hand-washing and surgical hand disinfection

Hand-washing according to EN 1500



2ND STEP Rub right palm over back of left hand and left palm over back of right hand



5TH STEP Rub the right and left thumbs



3RD STEP Rub palm against palm with interlocked, spread fingers



6TH STEP Rub closed tips of finger in the right and left palms

*Customizable software

The software that guides the process can be adapted to specific institutional regulations and conditions



SENSO Clean 1400

Elegant and ergonomic

The scrub-up trough is made of high-quality stainless steel and has been welded without seams in a special process The SENSO *Clean's* shape has been ergonomically designed for use and thus fulfills the requirements of the clinic staff and doctors exceptionally well

The SENSO Clean has no corners or edges, only curves. The risk of germs and bacteria collecting in poorly accessible points is greatly reduced

INTEGR©MED Integrated Medical Solutions

SENSO Clean 1400

Cross-functional intelligent synchronization



Room 1, 08:15

Room 2, 15:45



As an option, the system can transfer all data using a network Statistics on hand connection hygiene compliance can be compiled and divided by user groups (e.g. departments, hospitals, regions, etc.) **Operation Room** INTUITIVE Console ((•)) RFID Status Log in Hygiene control Protocol DICOM Server Audit Free Document added to patient files

INTEGR MED

SENSO Clean 1400

Mechanical UV/ozone disinfection ensures germ-free water quality



The backwash filter is maintenance-free and does not have to be changed^{*}

* The filter insert of backwash filters has to be replaced every 5 years.

The SENSO Clean UV/ozone filter system has been developed to establish a new way of chemical free and reliable disinfection - without changing the original water composition

Germs and microorganisms are neutralized by ozone and UV. The water is disinfected without leaving any residue

> A UV/0₃ disinfection of the process water reduces the presence of pathogenic germs in the water by 99.999%

The feed water is filtered by an integrated mechanical filter system with a maintenance-free backwash filter

INTEGROMED

INTEGROMED is your first point of contact when it comes to high-quality integrated solutions for the construction and furnishing of medical facilities. We capitalize on versatile, functional and at the same time attractive systems to create a unique environment for medical staff and patients. Not only do we deliver individual products, but also complete, tailor-made solutions.

INTEGROMED stands for competence, cooperation, responsibility and innovation.

The intelligent **SENSO** *Clean* is a genuine innovation in the field of medical hand-washing. The device has been developed on the basis of scientific research and the practical experience - for optimal hand hygiene in operating theaters.

The system features a new ergonomic design with automatic and manual operation modes. The completely touchless **SENSO** *Clean* can be used with a variety of filter systems, ranging from mechanical filters to sterile, antibacterial filters, depending on the clients' needs. The specially developed water quality control (WQC) system monitors the status of the filter system and informs the client of any maintenance requirements. Safety is guaranteed by a three-stage anti-scalding system. The large LCD front display shows vital information such as the steps for the hand-washing and disinfection process, time, water temperature, and the filling level of dispensers.

Furthermore, a UV quality control system to help train young doctors and verify that their hands have been properly decontaminated can also be integrated.



we believe in we

TECHNICAL DATA

Intended use	Scrub-up sink trough for hospitals with soap dispenser, disinfectant dispenser and UV module to ensure the personal hygiene of the hospital staff
Material	Glass, stainless steel
Display	10.1", 1,366 x 768 pixels, True Color, one for each washbasin
Display functions	Soap/disinfectant level, remaining time washing process, washing/disinfection movement, time, date, temperature, UV light
UV module	For an optical check of the disinfection success, 4 UV LEDs, wavelength 365nm, touchless switching
Water filter	Backwash filter
Hygiene rinsing	Adjustable hygiene program, after 12h of inactivity for soap, disinfectant and water
Pipe connection	3/4" hot and cold water
Temperature range	20.042.0°C, adjustable temperature limit
Operation	Touchless via IR sensors
Volume of soap and disinfectant tank	51
Lamp	LED 4 X 3.5 W
Power supply	230V/50Hz
Dimensions	Width: 700, 1,400, 2,100mm, Depth: 700 mm, Height: 1,500mm
Fastening	Wall mounting
Conformity	RoHS-compliant (lead-free), CE standard, 2014/30/EU EMC Directive, 2014/35/EU Low-Voltage Directive

SENSO Clean 700 / 1400 / 2100

